
Quarter 3 Update

Rohan Bhansali Avi Komarlingam

Data Processing











Bilinear Interpolation



Nearest Neighbor Interpolation



Bilinear Interpolation






(2539, 2705)

(512, 512)

Pixel Normalization






Min: 0.000, Max: 255.000


Min: 0.000, Max: 1.000

AP CNNs




DenseNet


DenseNet Fast DenseNet-121 DenseNet-161 DenseNet-169 DenseNet-201



DenseNet

- Alleviates the vanishing-gradient problem
- Strengthens feature propagation and reuse
- Substantially reduces the number of parameters





Model Summaries

DenseNet-121

DenseNet-161

Layers: 121 Layers: 161

Epochs: 3 Epochs: 3

Batch size: 8 Batch size: 5

Loss function: categorical cross entropy

Loss function: categorical cross entropy

Accuracy: 92.30% Accuracy: 92.15%

Time: 2321s/epoch Time: 3393s/epoch

"Results"

0	0.79	0	0	0	0	0	0	0	0	0.18	0.01	0.01
0	0.77	0	0	0	0	0	0.01	0	0	0.19	0.01	0.01
0	0.77	0	0	0	0	0	0.01	0	0	0.17	0.01	0.04
0	0.75	0	0	0	0	0	0.02	0	0	0.19	0.02	0.02
0	0.62	0	0	0	0	0	0	0	0	0.29	0	0.1
0	1	0	0	0	0	0	0	0	0	0	0	0
0	0.74	0	0	0	0	0	0	0	0	0.16	0.05	0.05
0	0.71	0	0	0	0	0	0.01	0	0	0.23	0.02	0.02
0	0.72	0	0	0	0	0	0.02	0	0	0.18	0.03	0.04
0	0.78	0	0	0	0	0	0	0	0	0.19	0.01	0.02
0	0.9	0	0	0	0	0	0	0	0	0.1	0	0
0	0.76	0	0	0	0	0	0.03	0	0	0.22	0	0
0	0.82	0	0	0	0	0	0.09	0	0	0.09	0	0

															1.0
Atelectasis -	0	0	0	0	0	0	0	0		0	0	0	0		
Cardiomegaly -	0	0	0	0	0	0	0	0		0	0	0	0		
Consolidation -	0	0	0	0	0	0	0	0		0	0	0	0		0.8
Edema -	0	0	0	0	0	0	0	0		0	0	0	0		
Enlarged Cardiomediastinum -	0	0	0	0	0	0	0	0		0	0	0	0		
Fracture -	0	0	0	0	0	0	0	0		0	0	0	0		0.6
- Lung Lesion -	0	0	0	0	0	0	0	0		0	0	0	0		
Lung Opacity -	0	0	0	0	0	0	0	0	1	0	0	0	0		
No Finding -		0	0	0	0	0	0	0	1	0	0	0	0		0.4
Pleural Effusion -	0	0	0	0	0	0	0	0	1	0	0	0	0		
				0	0	0	0	0	1	0	0	0	0		- 0.2
Pleural Other -	0	0	0												0.2
Pneumonia -	0	0	0	0	0	0	0	0		0	0	0	0		
Pneumothorax -	0	0	0	0	0	0	0	0		0	0	0	0		
	Atelectasis -	Cardiomegaly -	Consolidation -	Edema -	Enlarged Cardiomediastinum -	Fracture -	Lung Lesion -	Lung Opacity -	No Finding -	Pleural Effusion -	Pleural Other -	Pneumonia -	Pneumothorax -		- 0.0
Predicted label															


0	0.76	0	0	0.02	0.04	0.03	0.03	0.01	0	0.02	0.06	0.02
0	0.78	0	0	0.03	0.03	0.04	0.02	0.02	0	0.03	0.05	0.01
0	0.74	0	0	0.03	0.06	0.06	0.01	0.02	0	0.03	0.05	0.02
0	0.8	0	0.01	0.01	0.02	0.02	0.04	0.01	0.01	0.02	0.05	0.02
0	0.71	0	0	0	0	0.05	0	0.05	0	0	0.19	0
0	0.5	0	0	0	0.5	0	0	0	0	0	0	0
0	0.95	0	0	0	0	0	0	0	0	0	0.05	0
0	0.77	0	0	0.02	0.03	0.01	0	0.01	0	0.04	0.07	0.03
0	0.71	0	0	0.01	0.05	0.03	0.02	0.02	0	0.02	0.12	0.02
0	0.76	0	0	0.03	0.03	0.01	0.04	0.01	0.01	0.05	0.06	0.01
0	0.7	0	0	0	0	0	0	0	0	0.1	0.1	0.1
0	0.81	0	0	0.03	0	0	0.03	0.03	0	0	0.08	0.03
0	0.73	0	0	0	0	0	0.09	0.09	0	0	0	0.09

Possible Causes

1. Domination of one class over the others due to

overrepresentation

Equalizing Class Sizes




Possible Causes

 Domination of one class over the others due to overrepresentation

2. Errors in preprocessing images

Pixel Denormalization




Min: 0.000, Max: 1.000

Min: 0.000, Max: 255.000

Possible Causes

- Domination of one class over the others due to overrepresentation
- 2. Errors in preprocessing images
- Issues with model architecture

Different Model Architectures



Even More Different Model Architectures

kNN

Logistic Regression

SVM

	Predicted											
		Enlarged Cardiomediastinum	Lung Lesion	Lung Opacity	Σ							
	Enlarged Cardiomediastinum	51.7 %	38.7 %	9.6 %	551							
Actual	Lung Lesion	36.0 %	55.3 %	8.8 %	684							
Act	Lung Opacity	46.9 %	43.0 %	10.1 %	435							
	Σ	735	778	157	1670							
		Predicted										
		Enlarged Cardiomediastinum	Lung Lesion	Lung Opacity	Σ							
	Enlarged Cardiomediastinum	43.0 %	35.6 %	21.4 %	551							
nal	Lung Lesion	27.2 %	52.5 %	20.3 %	684							
Actual	Lung Opacity	32.9 %	38.4 %	28.7 %	435							
	Σ	566	722	382	1670							
		Predicted										
		Enlarged Cardiomediastinum	Lung Lesion	Lung Opacity	Σ							
	Enlarged Cardiomediastinum	53.9 %	22.7 %	23.4 %	551							
nal	Lung Lesion	51.3 %	29.2 %	19.4 %	684							
Actual	Lung Opacity	54.5 %	19.5 %	26.0 %	435							
	Σ	885	410	375	1670							

Scientist Outreach





Future Work

Train and optimize AP, PA, and lateral CNNs

Fuse networks via DualNet architecture

Incorporate patient information (syndromes, history, etc.) into diagnosis